

40G fiber optic data interface

4 images in 10 ns 8 images in 580 ns **exposure time 2.5 ns** with 18 mm intensifier









As an expansion of our well established pco.dicam C4 family, we are now adding the pco.dicam C4 UHS models, which are based on ultra fast 18 mm image intensifiers. Similar to their 25 mm counterparts they are available with a range of high quality photo cathodes matching different application requirements. The pco.dicam C4 UHS is the first multi-channel intensified camera system which exploits the full performance inherent to **scientific CMOS** sensor technology.

With its high-end optical beam splitters you are able to equally distribute the input light to the 4 image intensifiers. They are coupled with the pco.dicam C1 proven tandem lenses to the 16 bit 2.3 MPixel sCMOS sensor. It is the most flexible configuration of 8 individual exposure times and their corresponding interframing times, which makes the camera so unique. The 40G fiber optic based data interface (CLHS FOL) guarantees you uncompressed and robust 16 bit data transfer of 572 full frames per second via optical fiber over virtually any distance.

features&benefits

572 fps @ full 2.3 MPixel resolution	high frame rates at high resolution for imaging of dynamic processes
1.1 e- readout noise	lowest readout noise of any gated intensified camera system
16 bit digitization	taking advantage of the higher dynamic range possible from high end image intensifiers
optical coupling via ultra-speed tandem lens	outstanding image quality with high transmission efficiency and no artifacts
tandem lens with 0.53 : 1 image scaling	full 18 mm diameter of intensifier output is imaged (lossless) onto an sCMOS sensor
40G fiber optic based data interface	fiber optic interface virtually covers any distance without deploying additional interface converters or signal amplifiers with immunity to EMI
4 x 660 MByte/s image data rate	highest sustained image data rate of any intensified camera system on the market; no limitations for recording duration
double shutter mode with 300 ns interframing time	two consecutive full resolution images with a configurable minimum interframing time of 300 ns on each of the 4 channels
2.3 MPixel sCMOS sensor	overcomes CCD limitations in terms of speed and sensitivity
enhanced extinction ratio gating	fast MCP gating for improved extinction ratio for the blue and uv part of the spectrum
additional optical trigger input	robust trigger transmission over long distance in EMC critical environments
lense remote controller (optional)	convenient remote lens control for camera systems inaccessible during an experiment
selected highly homogeneous image intensifiers	uses best image intensifier quality available on the market
50 ns trigger to exposure start delay	ultra-fast camera reaction to trigger event
2.5 ns gating with 18 mm intensifier	captures fast transient phenomena
extensive and highly precise IN/ OUT signaling	allows for perfect synchronization in any experimental setup as timing master or slave
configurable delay in steps of 1 ns	flexible adaptation to synchronization needs

945 East 11th Avenue Tampa, FL 33605

Phone: (813) 984-0125

Contact: Sales@ pyramidimaging.com

https://pyramidimaging.com




camera components overview

>> top view

front part with optical input and beam splitter optics

- A collimator lens generates bundles of parallel rays with focus set to infinity.
- In between, 3 double prisms provide a 50:50 beam distribution under a 90° angle.
- 3 2 single prisms act as 99.9 % reflection mirrors.
- Spectral filters can be mounted individually for each of the 4 light channels (CH1, CH2, CH3, CH4).
- The imaging lenses of each channel focus the parallel bundles onto the photocathode of the image intensifier.
- 6 Image intensifier (See more information on the next page.)





945 East 11th Avenue Tampa, FL 33605


Phone: (813) 984-0125

Contact: Sales@ pyramidimaging.com



camera components overview

- 6 image intensifier
- 7 optical coupling lens system
- 8 sCMOS image sensor
- 9 camera system
- 10 10G fiber optic based interface








technical specifications

>> image intensifier (4x)

type	HighRes micro channel plate (MCP) 6 μm channel
input window	synthetic silica, borosillicate
photocathode material	S20, GaAs, GaAsP (others on request)
image intensifier pitch distance	6 μm
image intensifier MCP type	single stage low resistance MCP for high strip current
MCP operational modes	continuous gated for enhanced extinction ratio
image intensifier diameter	18 mm
phosphor screen material	P43, P46
output window	glass
image intensifier system resolution	> 50 lp/mm @ 5 % MTF typical (depends on phosphor)
shortest gating time	2.5 ns

» image intensifier input window

Typical transmittance of image intensifier input window materials.



data courtesy of Hamamatsu Photonics

Due to the optical properties of the beam-splitter optics, there is no uv transmission below 380 nm. Intensifiers with MgF2 input window are not available. Standard input window for S20 photocathodes is synthetic silica.

GaAs and GaAsP photocathodes are deposited on borosilicate glass.




945 East 11th Avenue Tampa, FL 33605



>> image intensifier photocathode characteristics

Spectral sensitivities of different photocathode materials: S20 (multialkali), GaAs, GaAsP



data courtesy of Hamamatsu Photonics

photocathode material	peak wavelength [nm]	typical quantum efficiency at peak wavelength [%]	dark counts [s ⁻¹ /cm ²]
S20 (multialkali)	250	20	1500
GaAs	650	30	30,000
GaAsP	500	55	10,000

data courtesy of Hamamatsu Photonics

>> image intensifier phosphor

phoophor	phosphor decay (typ.) to		peak	typical
phosphor	10 %	1 %	emission	efficiency
P43	1 ms	4 ms	545 nm	100 %
P46	0.2 - 0.4 µs	2 µs	530 nm	30 %

You can combine all photocathode materials with P43 or P46 phosphor. Whereas the P43 phosphor has a much brighter emission than the P46 phosphor, it has a rather long decay time, i.e. the time required until the phosphor emission fades out after the excitation by electron bombardement has been stopped. This decay time is therefore critical for fast image repetition rates primarily in double image application or when operating the camera in spectroscopic mode with line rates in the kHz range.








technical specifications

>> optical coupling lens system of the detector units (4x)

ultra-speed tandem lens between image intensifier & sCMOS

transmission efficiency	> 30 %
vignetting	< 3 %
resolution	> 60 lp/mm
scaling rates	B=0.53 for 18 mm intensifier



The projected image circle is completely covered by $1504 \times 1504 \times 1504$





technical specifications

>> sCMOS image sensor

Each detector unit of this unique 4 channel design is equipped with an sCMOS image

type of sensor	scientific CMOS (sCMOS)
resolution (h x v)	1504 x 1504 active pixel
pixel size (h x v)	6.5 µm x 6.5 µm
sensor format / diagonal	9.8 mm x 9.8 mm / 13.8 mm
shutter mode	single image double image
MTF ¹	76.9 lp/mm (theoretical)
fullwell capacity	15,000 e- for P46 phosphor 30,000 e- for P43 phosphor
readout noise ²	1.1 med / 1.5 ms e ⁻ single image 2.2 med / 2.5 ms e ⁻ double image
dynamic range	13,600 : 1 (82.7 dB) for P46 phosphor 27,200 : 1 (88.7 dB) for P43 phosphor
quantum efficiency	58 % for P43 peak emission @ 545 nm 57 % for P46 peak emission @ 530 nm
spectral range	300 nm 1000 nm
dark current ³	< 0.6 e ⁻ /pixel/s @ 7 °C
DSNU	1.0 e ⁻ ms
PRNU	< 0.6 %
anti blooming factor	1:10,000

>> frame rate table4

	C1 UHS	C4 UHS	C8 UHS
1504 x 1504	143 fps	572 fps	1144 fps
1504 x 1024	210 fps	840 fps	1680 fps
1504 x 512	414 fps	1656 fps	3312 fps
1504 x 256	807 fps	3228 fps	6456 fps
1504 x 128	1535 fps	6140 fps	12,280 fps
1504 x 64	2795 fps	11,180 fps	22,360 fps
1504 x 32	4739 fps	18,956 fps	37,912 fps
1504 x 16	7266 fps	29,064 fps	58,128 fps
1280 x 1024	210 fps	840 fps	1680 fps
640 x 480	441 fps	1764 fps	3528 fps
320 x 240	858 fps	3432 fps	6864 fps

Pyramid Imaging



¹ Modulation transfer function.
2 The readout noise values are given as median (med) and root mean square (rms) values due to the different noise models, which can be used for evaluation.

All values are raw data without any filtering. Measurements with dark current compensation. Exposure time < 1 μ s.




frame rates

Due to the special 4 channel design of the pco.dicam C4 and the flexible timing possibilities, extremely high burst frame repetition rates are feasible. In single image mode you can record sequences of 4 ultra fast images and in double image mode sequences of 8 ultra fast images. Examples for such extreme frame repetition rates are given

>> continuous imaging

2048 x 2048	572 fps
2048 x 256	up to 3200 fps
2048 x 16	> 27,000 fps
single image mode	4 images of 2.5 ns exposure time with 0 ns interframing time: 400,000,000 fps This 4 image sequence can be repeated every 7 ms
double image mode	8 images of 20 ns exposure time with 60 ns interframing time: 12,500,000 fps This 8 image sequence can be repeated every 28 ms

>> perfect fit: phosphor emission vs. sCMOS quantum efficiency



This chart describes the spectral situation for the internal imaging of the image intensifier's phosphor output screen to the sCMOS sensor of the camera detector module. This imaging is done by the highly efficient tandem lens system.

Please note: The spectral sensitivity relevant for your experiment is solely determined by the QE curve of the photocathode material of the image intensifier (page 6).



945 East 11th Avenue Tampa, FL 33605





technical specifications

>> detector unit (4x)

frame rate	143 fps @ 1504 x 1504 pixel > 7000 fps @ 1504 x 16 pixel
dynamic range A/D ⁵	16 bit
pixel scan rate	286.0 MHz
binning horizontal	x1, x2, x4
binning vertical	x1, x2, x4
region of interest (ROI)	horizontal: steps of 4 pixels vertical: steps of 1 pixel
non linearity	<1%
cooling method	+ 7 °C stabilized, 1 stage peltier with forced air (fan)
input signals	electrical trigger, arm input (TTL level, BNC connectors), gate disable (high-speed TTL input, BNC connectors)
output signals	gate/expos out monitor, user monitor output (TTL level, BNC connectors)
time stamp	in image (1 μs resolution)

>> exposure modes

single image mode

exposure times	fixed values < 20 ns depending on model (e.g. 2.5, 4,, 10 ns), 20 ns 250 ns (1 ns steps), 250 ns 1 s (10 ns steps)
delay times	0 ns 250 ns (1 ns steps), 250 ns 1 s (10 ns steps)
maximum repetition frequency	200 kHz sustained, 3.3 MHz burst
insertion delay	
trigger input to exposure out	19 ns
trigger input to optical open	49 ns
jitter	
trigger input to exposure out	35 ps rms
trigger input to optical open	150 ps rms
double image mode	
exposure times	20 ns 1 ms (in 10 ns steps)
delay settings	0 ns 10 ms (in 10 ns steps)
interframing time on every channel	300 ns 10 ms (in 10 ns steps)

⁵ The high dynamic signal is simultaneously converted at high and low gain by two 11 bit A/D converters and the two 11 bit values are sophistically merged into one 16 bit value.



945 East 11th Avenue Tampa, FL 33605



» general camera system

power supply	110 - 230 V
power consumption	180 W
weight	43.3 kg
operating temperature	+ 10 °C + 40 °C
operating humidity range	10 % 80 % (non-condensing)
storage temperature range	- 10 °C + 60 °C
optical mount	F-mount optional: Canon EF mount
lens remote controller (optional)	electronic control for Canon EF lenses
maximum cable length	10 km (CLHS FOL)
input signals	master trigger electrical and optical
CE / FCC certified	yes





technical specifications

>> camera interface (4x)

data transfer	Camera Link HS, FOL (Single F2, 1X1, S10) single 4 port frame grabber for PCI Express
maximum cable length	10 km (CLHS FOL)
master input signals	optical trigger (FOL), electrical trigger, arm input (TTL level, BNC connectors)
additional input signals per channel	electrical trigger, arm input (TTL level, BNC connectors), gate disable (high-speed TTL input, BNC connectors)
additional output signals per channel	gate/expos out monitor, user monitor output (TTL level, BNC connectors)





945 East 11th Avenue Tampa, FL 33605




technical specifications

>> lens remote controller

The optional Canon lens controller enables you to connect electronic EF and EF-S Canon lenses allowing to remote control focus and aperture of those lenses.

» dimensions

Camera equipped with F-mount lens adapter. All dimensions are given in millimeter.



>> camera view





945 East 11th Avenue Tampa, FL 33605

Phone: (813) 984-0125

Contact: Sales@ pyramidimaging.com

https://pyramidimaging.com



>> applications

laser induced incandescence (LII) | shock wave physics | laser induced breakdown spectroscopy (LIBS) particle image velocimetry (PIV) | time resolved spectroscopy | plasmaphysics | laser induced fluorescence (LIF) ballistics | combustion

≫ software



With pco.camware you control all camera settings, the image acquisition, and the storage of your image data. The pco.sdk is the complementary software development kit. It includes dynamic link libraries for user customization and integration on Windows PC platforms. Drivers for popular third party software packages are also available for

All these items like pco.camware, pco.sdk, and third party drivers are free to download at www.pco.de

third party integrations











945 East 11th Avenue Tampa, FL 33605

Phone: (813) 984-0125

Contact: Sales@ pyramidimaging.com



customization

» possible combinations

photocathode	input window	phosphor
S20 selected	synthetic silica	P46
		P43 ⁶
GaAs standard	borosilicate	P46
		P43 ⁶
GaAsP standard	borosilicate	P46
		P43 ⁶

 $^{^{\}mbox{6}}$ P43 phosohor cannot be used, if 8 fast images with interframing times < 1 ms are required.

18 mm image intensifiers are available in two quality grades.

Standard	quality specified for central 13.5 mm x 10 mm square region corresponding to 1100 x 810 pixel sCMOS sensor resolution
Selected	quality specified for 18 mm diameter area corresponding to full 1504 x 1504 pixel sCMOS sensor resolution, extinction ratio 10 times higher than standard grade, image intensifiers with S20 photocathode exclusively come in selected grade quality, contact our technical sales team for further details on the two quality grades

>> select optical mount

F-mount

Canon EF mount

>> select interface

type of fiber optic interface (CLHS FOL) module in camera and frame grabber

SM SFP+ up to 10 km $\,$

MM SFP+ up to 300 m

FOL cable length default: 10 m



945 East 11th Avenue Tampa, FL 33605

Phone: (813) 984-0125

Contact: Sales@ pyramidimaging.com

https://pyramidimaging.com

contact

pco europe +49 9441 2005 50 info@pco.de pco.de

pco america

+1 866 678 4566 info@pco-tech.com pco-tech.com

pco asia +65 6549 7054 info@pco-imaging.com pco-imaging.com

pco china

+86 512 67634643 info@pco.cn pco.cn



for application stories please visit our website



An Excelitas Technologies Brand